

**PHILIPS**

Ultrasound

# Noninvasive liver fibrosis assessment

ElastPQ ultrasound  
shear wave elastography

Richard G. Barr, MD, PhD, FACR, Diagnostic Radiology, Southwoods Imaging, Youngstown, Ohio.

## Liver disease – cause and prevalence

Diffuse liver disease is one of the major health problems in the world. It can result from many causes, including viral hepatitis (hepatitis B or hepatitis C), non-alcoholic or alcoholic fatty liver disease, autoimmune hepatitis, drug-induced liver injury, primary biliary cirrhosis, and several other less frequent etiologies. It is estimated that 360 million and 180 million people worldwide are infected with viral hepatitis B and C respectively. Between 500,000 and 700,000 people die annually as a result of hepatitis B virus infection, and more than 350,000 people are estimated to die each year from hepatitis C-related liver disease.<sup>1-3</sup> Chronic liver damage results in hepatic fibrosis, characterized by an increase in extracellular matrix material produced by fibroblast-like cells. This process results in liver fibrosis that can progress to cirrhosis with distortion of normal liver architecture and portal hypertension. According to SRU consensus, in select patients, elastography may eliminate the need for liver biopsy for staging fibrosis.<sup>4</sup>



ElastPQ is available on  
Philips EPIQ 7, EPIQ 5,  
Affiniti 70 and iU22.

# Noninvasive liver fibrosis assessment

## ElastPQ ultrasound shear wave elastography

### Diagnosis and staging of liver disease

Accurately staging the degree of liver fibrosis is extremely important to determine if antiviral therapy is appropriate, and to predict treatment outcome and malignant potential. With current drug therapy, early stage fibrosis may be reversible.<sup>5</sup>

The histologic evaluation of liver biopsies is carried out using scoring systems that produce values for various categories of inflammation (grade), and fibrosis (stage). There are several scoring systems, all categorizing similar features. In the assessment of chronic HCV hepatitis, the most reproducible scoring system is the METAVIR. On the METAVIR scoring system, liver fibrosis is evaluated semi-quantitatively and staged on a five-point scale from 0 to 4 (F0: absent; F1: enlarged fibrotic portal tract; F2: peri-portal or initial portal-portal septa but intact architecture; F3: architectural distortion but no obvious cirrhosis; and F4: cirrhosis).<sup>6</sup>

The gold standard for diagnosis and staging of liver fibrosis has been liver biopsy. In addition to being an invasive procedure with potential complications of bleeding and severe pain, sampling error is an intrinsic problem due to the small sample size in a heterogeneous process.<sup>7,8</sup> Inter-observer variability also limits diagnostic consistency.<sup>9-11</sup> The development of several blood markers such as platelets, hyaluronic acid, type IV collagen, aminotransferase/platelet ratio index (APRI) and algorithm based serum models (Fibro Index, FIB-4, and Fibro Test) have been used but are affected by factors unrelated to the liver.

ElastPQ is available on Philips EPIQ 7, EPIQ 5, Affiniti 70 and iU22.

### ElastPQ – a new era in liver disease assessment

A technique called ElastPQ uses ultrasound shear wave elastography (ARFI technology) to provide a noninvasive, reproducible, and easily performed method of assessing liver fibrosis. A special pulse sequence technique that uses existing transducers produces shear waves in tissue and then measures the propagation speed of the waves. Now liver stiffness samples can be acquired during a routine ultrasound examination of the liver. According to a recent study, using shear wave elastography may help reduce or avoid conventional liver biopsies.<sup>12</sup> Instead of a costly and painful biopsy procedure, an easy ultrasound exam becomes the routine method to assess liver disease status.

### Suggested thresholds using ElastPQ shear wave elastography in patients with hepatitis C<sup>4</sup>

Liver fibrosis staging	kPa	m/s
No clinical significant fibrosis: METAVIR stage equal or less F2 Unlikely to need follow-up	<5.7	<1.37
Advanced fibrosis and/or cirrhosis: METAVIR stage of F4 and some stages of F3 – clinically significant fibrosis	>15	>2.2



## Performing an ElastPQ shear wave examination

With increasing fibrosis, the liver becomes stiffer, which can be monitored using shear wave elastography.<sup>13,14</sup> With this technique, during an ultrasound exam after selecting the best acoustic window, a region of interest (ROI) is placed in an area of the liver perpendicular to the liver capsule, taking care not to include large vasculature or biliary structures. An intercostal imaging approach targeting segments 7 or 8 of the liver has been shown to provide more reliable measurements. Serial measurements are taken while the patient suspends respiration and a report is generated. Literature suggests that minimum of 10 measurements should be obtained out of which 60% should be “good”.<sup>4</sup> The average of these measurements is then used to estimate the degree of liver stiffness and correlate with a predicted biopsy METAVIR score.

Although this technique typically shows strong correlation, there are several confounding factors that may distort results, such as liver inflammation, liver congestion, and biliary obstruction. In some cases, distinguishing normal from very mild disease may be difficult; also, some moderate and severe disease may look similar.

ElastPQ uses ultrasound shear wave elastography to provide a **noninvasive, reproducible, and easily performed** method of assessing liver fibrosis.

## The potential for change in clinical practice and reduction in medical costs

The United Kingdom’s NICE (National Institute for Health and Care Excellence) updated their clinical guidelines for managing patients with chronic viral hepatitis B infections.<sup>15</sup> The newest guidelines recommend that transient elastography is offered as the initial test for liver disease in adults newly referred for assessment, and annual reassessment of liver disease using transient elastography is offered to adults who are not taking antiviral treatment. The guidelines only target liver biopsy in very specific conditions. This provides a less invasive approach to managing liver disease and offers the opportunity to control medical costs.

The European Federation of Societies for Ultrasound in Medicine and Biology and World Federation for Ultrasound in Medicine and Biology have produced guidelines for the use of elastography for evaluation of liver stiffness.

When compared to transient elastography exams, Philips ElastPQ shear wave elastography has been found to show similar accuracy and is an image-guided exam that shows precisely where stiffness measurements are taken.<sup>4,12</sup> ElastPQ has also been shown to allow measurements in patients with high BMI (Body Mass Index) and patients with ascites.

As clinical practice standards of care adapt to appropriate use guidelines, ElastPQ elastography may offer an ideal way to routinely monitor liver tissue stiffness, and reduce other, more costly and invasive methods of testing. CMS CPT add-on codes<sup>16</sup> for provider tracking of elastography exams exist currently and provide the preliminary steps for future reimbursement.

## Key take-aways for ElastPQ liver stiffness assessment

- Easily combine a routine ultrasound imaging exam of the liver anatomy with targeted tissue stiffness values
- Assess liver fibrosis in patients with clinically suspected disease even before abnormalities are detected with ultrasound imaging
- Evaluate and obtain a baseline stiffness value in patients with chronic liver disease
- Follow up patients under treatment to monitor progression, stabilization or regression of liver disease
- Help avoid the need for liver biopsies when elastography results are consistent with other clinical findings

## References

- 1 European Association for the Study of the L. EASL Clinical Practice Guidelines: management of hepatitis C virus infection. *Journal of Hepatology*. 2011;55(2):245-64.
- 2 European Association for the Study of the L. EASL Clinical Practice Guidelines: management of chronic hepatitis B virus infection. *Journal of Hepatology*. 2012;57(1):167-85.
- 3 World Health Organization. Viral hepatitis. Report from the Secretariat. Sixty-third World Health Assembly. 2010.
- 4 Barr RG, Ferraioli G, Palmeri ML, Goodman ZD, Garcia-Tsao G, Rubin J, Garra B, Myers RP, Wilson SR, Rubens D, Levine D. Elastography Assessment of Liver Fibrosis: Society of Radiologists in Ultrasound Consensus Conference Statement. *Radiology*, Sep 2015, Vol. 276: 845-861, 10.1148/radiol.2015150619
- 5 Sohrabpour AA, Mohamadnejad M, Malekzadeh R. The Reversibility of Cirrhosis. *Disclosures Aliment Pharmacol Ther*. 2012;36(9):824-832.
- 6 Bedossa P, Poynard T. An algorithm for the grading of activity in chronic hepatitis C. The METAVIR Cooperative Study Group. *Hepatology*. 1996;24(2):289-93.
- 7 Bravo AA, Sheth SG, Chopra S. Liver biopsy. *The New England Journal of Medicine*. 2001;344(7):495-500.
- 8 Cadranet JF, Rufat P, Degos F. Practices of liver biopsy in France: results of a prospective nationwide survey. For the Group of Epidemiology of the French Association for the Study of the Liver (AFEFL). *Hepatology*. 2000;32(3):477-81.
- 9 Maharaj B, Maharaj RJ, Leary WP, et al. Sampling variability and its influence on the diagnostic yield of percutaneous needle biopsy of the liver. *Lancet*. 1986;1(848 0):523-5.
- 10 Bedossa P, Dargere D, Paradis V. Sampling variability of liver fibrosis in chronic hepatitis C. *Hepatology*. 2003;38(6):1449-57.
- 11 Regev A, Berho M, Jeffers LJ, et al. Sampling error and intraobserver variation in liver biopsy in patients with chronic HCV infection. *The American Journal of Gastroenterology*. 2002;97(10):2614-8.
- 12 Ferraioli G, Tinelli C, Lissandrin R, et al. Reproducibility and performance of a new point shear wave elastography technique for assessing fibrosis in chronic hepatitis C. *World Journal of Gastroenterology*. 2014;20(16):4787-4796.
- 13 Ferraioli G, Lissandrin R, Zicchetti M, Filice C, Calliada F CM, Ferraiolo G. Sono-Elastography: Main Clinical Applications in Diffuse Liver Diseases. Edimes Srl Edizioni Medico Scientifiche – Pavia, ISBN 978-88-88541-09-9; March 30, 2013.
- 14 Yu H, Wilson SR. New noninvasive ultrasound techniques: can they predict liver cirrhosis? *Ultrasound Quarterly*. 2012;28(1):5-11.
- 15 NICE clinical guideline 165, guidance.nice.org.uk/cg165.
- 16 AMA Website: <http://www.amaassn.org/resources/doc/cpt/cptcat3codes.pdf> for implementation in 2014.

